

[Cedrus](#)[Products](#)[Ordering](#)[Support](#)[XID Home](#)[Device Properties](#)[Timing Features](#)[Supported Protocols](#)[Accessory Connector](#)[Voice Key](#)[Test Procedures](#)[Misc. Commands](#)

XID Accessory Connector Command

Set Output Logic Mode (a0)

What it Does

Sets the MRI scanner's trigger debouncing time.

Using This Command

Send bytes **f5** followed by a single byte indicating the debounce time in milliseconds.

Corresponding Inquiry Command

Send bytes **_f5**. The XID device returns **_f5** followed by a byte indicating the debounce time.

© Copyright 2014 **Cedrus Corporation**, P.O. Box 6309, San Pedro, CA 90734 - USA

Toll Free: (800) 233-7871, other: (310) 548-9595. Send us [an email](#) or join [mailing list](#). 