

[Cedrus](#)[Products](#)[Ordering](#)[Support](#)[XID Home](#)[Device Properties](#)[Timing Features](#)[Supported Protocols](#)[Accessory Connector](#)[Voice Key](#)[Test Procedures](#)[Misc. Commands](#)

XID Accessory Connector Command

Set 'OR' Line (a0)

What it Does

Sets one of the Accessory Connector's I/O lines as a logical OR line.

On **Lumina LP-400** response pads: line AC6 is set as a logical "OR" line. It is turned on when any key is pressed on the response pads.

On **SV-1** voice key devices: this command has no effect.

Using This Command

To Turn Logical 'OR' Mode	Send Bytes
OFF	ao0
ON	ao1

Corresponding Inquiry Command

Send bytes **_ao**. The XID device returns **_ao0** or **_ao1**.