

[XID Home](#)
[Device Properties](#)
[Timing Features](#)
[Supported Protocols](#)
[Accessory Connector](#) ▶

[Voice Key](#)
[Test Procedures](#)
[Misc. Commands](#) ▶

XID Miscellaneous Commands

The following table provides a summary of the miscellaneous XID commands available:

Command Name	Command	Description
Set Serial Port Speed	f1	Sets the serial port speed
Set Locking Level	f2	Locks or unlocks the XID device's manual dials and other controls
Re-program Flash	f3	Makes the firmware ready to download a new version
Set Scanner Trigger Default	f4	Sets the MRI scanner's default value
Set Trigger Debounce Time	f5	Sets the MRI scanner's debounce time
Set Button Debounce Time	f6	Sets the push buttons' debounce time
Restore Factory Defaults	f7	Resets the XID device to factory defaults
Get Battery Level	_f8	Returns the battery level on battery-operated devices